

Precolumn types and replacement guidelines, etc.

List of precolumns for each analysis

Separation column Type Size	Sample type	Mobile phase pH	Target substance	Recommended precolumn		
				Types of fillers	Size	Position
CAX φ 2.0 x 200 mm	microdialysis	6	NE, DA, 5-HT	CA-ODS	φ 4 x 5 mm	Upstream from injector
PP-ODSIII φ 4.6 x 30 mm	microdialysis	6	DA, 5-HT	CA-ODS	φ 4 x 5 mm	Upstream from injector
SC-5ODS φ 2.1 x 150 mm	microdialysis	2~4	Monoamine and Metabolites, etc.	AC-ODS	φ 4 x 5 mm	Upstream from injector
	Other than microdialysis	2~4		AC-ODS	φ 3 x 4 mm	Just before the separation column
	Other than microdialysis	4~		CA-ODS	φ 3 x 4 mm	Just before the separation column
SC-5ODS φ 3.0 x 150 mm	Other than microdialysis	2~4	Monoamine and Metabolites, etc.	AC-ODS	φ 4 x 5 mm	Just before the separation column
	Other than microdialysis	4~		CA-ODS	φ 4 x 5 mm	Just before the separation column
CA-5ODS φ 2.1 x 150 mm	microdialysis	5~6	Monoamine	CA-ODS	φ 4 x 5 mm	Upstream from injector
	Other than microdialysis	5~6		CA-ODS	φ 3 x 4 mm	Just before the separation column
AC-GEL φ 2.0 x 150 mm	microdialysis		Acetylcholine	CH-GEL	φ 4 x 5 mm	Upstream from injector
	Other than microdialysis			CH-GEL	φ 3 x 4 mm	Just before the separation column
AC-GEL φ 4.6 x 150 mm	Other than microdialysis		Acetylcholine	CH-GEL	φ 4 x 5 mm	Just before the separation column
GU-GEL φ 4.6 x 150 mm	microdialysis		Glutamate	CH-GEL	φ 4 x 5 mm	Upstream from injector
	Other than microdialysis			CH-GEL	φ 4 x 5 mm	Just before the separation column

Guidelines for precolumn replacement

1. Microdialysis sample → Approx. **500** samples
2. Tissue homogenized sample for monoamine analysis → Approx. **100 - 200** samples
3. Tissue homogenized sample for acetylcholine analysis → Approx. **50 - 100** samples
4. Blood sample → Approx. **100 - 200** samples

In other cases, refer to items 1 to 4 above.

Also, if the pressure exceeds **1MPa** in the precolumn only, it is recommended to replace it.

List of precolumn models

Types of fillers	φ 3 x 4 mm	φ 4 x 5 mm
AC-ODS	PC-03-AC	PC-04-AC
CA-ODS	PC-03-CA	PC-04-CA
CH-GEL	PC-03-CH	PC-04-CH

